

## NFRKB Polyclonal Antibody

Catalog No: #27297

Package Size: #27297-1 50ul #27297-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	NFRKB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human NFRKB (NP_001137307.1).
Other Names	NFRKB; INO80G; nuclear factor related to kappaB binding protein
Accession No.	Swiss-Prot#:Q6P4R8NCBI Gene ID:4798
Uniprot	Q6P4R8
GeneID	4798;
Calculated MW	185kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using NFRKB at 1:1000 dilution.

## Background

---

Note: This product is for in vitro research use only