

GRHL2 Polyclonal Antibody

Catalog No: #27312



Package Size: #27312-1 50ul #27312-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GRHL2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human GRHL2 (NP_079191.2).
Other Names	GRHL2;BOM;DFNA28;ECTDS;TFCP2L3
Accession No.	Uniprot:Q6ISB3GenelD:79977
Uniprot	Q6ISB3
GenelD	79977
Calculated MW	75kDa
SDS-PAGE MW	71kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using GRHL2 antibody.



Western blot analysis of extracts of Rat lung, using GRHL2 antibody.

Background

The protein encoded by this gene is a transcription factor that can act as a homodimer or as a heterodimer with either GRHL1 or GRHL3. Defects in this gene are a cause of non-syndromic sensorineural deafness autosomal dominant type 28 (DFNA28).

Note: This product is for in vitro research use only