

CHD8 Antibody

Catalog No: #45883

Package Size: #45883-1 50ul #45883-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CHD8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	CHD8 Antibody detects endogenous levels of total CHD8
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CHD8
Target Name	CHD8
Other Names	CHD8, ATP-dependent helicase CHD8, Axis duplication inhibitor, CHD-8, HELSNF1, KIAA1564, Duplin, Helicase with SNF2 domain 1
Accession No.	Swiss-Prot#: Q9HCK8NCBI Gene ID:57680
Uniprot	Q9HCK8
GeneID	57680;
Calculated MW	290kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of lysates of HeLa cell and mouse brain, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only