

coiled-coil domain containing 102B antibody

Catalog No: #23151

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|-----------------------|--|
| Product Name | coiled-coil domain containing 102B antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Purified by antigen-affinity chromatography. |
| Applications | WB |
| Species Reactivity | Hu |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 244 and 509 of Human CCDC102B |
| Target Name | coiled-coil domain containing 102B |
| Accession No. | Swiss-Prot:Q68D86Gene ID:79839 |
| Uniprot | Q68D86 |
| GeneID | 79839; |
| Concentration | 1mg/ml |
| Formulation | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative. |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |

Application Details

Predicted MW: 60kd

Western blotting: 1:500-1:3000

Images



Sample (30ug whole cell lysate)
A: Hep G2
B: MOLT4
7.5% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only