

coiled-coil domain containing 102B antibody

Catalog No: #23151

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	coiled-coil domain containing 102B antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 244 and 509 of Human CCDC102B
Target Name	coiled-coil domain containing 102B
Accession No.	Swiss-Prot:Q68D86Gene ID:79839
Uniprot	Q68D86
GeneID	79839;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 60kd

Western blotting: 1:500-1:3000

Images



Sample (30ug whole cell lysate)
 A: Hep G2
 B: MOLT4
 7.5% SDS PAGE
 Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only