

COG6 Polyclonal Antibody

Catalog No: #27369

Package Size: #27369-1 50ul #27369-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	COG6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human COG6 (NP_065802.1).
Other Names	COG6; CDG2L; COD2; SHNS; component of oligomeric golgi complex 6
Accession No.	Swiss-Prot#:Q9Y2V7NCBI Gene ID:57511
Uniprot	Q9Y2V7
GeneID	57511;
Calculated MW	73kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:1000 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using COG6 at 1:1000 dilution.

Background

This gene encodes a subunit of the conserved oligomeric Golgi complex that is required for maintaining normal structure and activity of the Golgi apparatus. The encoded protein is organized with conserved oligomeric Golgi complex components 5, 7 and 8 into a sub-complex referred to as lobe B. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only