

SPSB1 Polyclonal Antibody

Catalog No: #27375

Package Size: #27375-1 50ul #27375-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SPSB1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human SPSB1
Other Names	SPSB1; SSB-1; SSB1; SPRY domain-containing SOCS box protein 1
Accession No.	Swiss-Prot#:Q96BD6NCBI Gene ID:80176
Uniprot	Q96BD6
GeneID	80176;
Calculated MW	31kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse heart, using SPSB1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only