

## FBXO38 Polyclonal Antibody

Catalog No: #27377

Package Size: #27377-1 50ul #27377-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FBXO38 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of human FBXO38
Other Names	FBXO38; Fbx38; HMN2D; MOKA; SP329; F-box protein 38
Accession No.	Swiss-Prot#:Q6PIJ6NCBI Gene ID:81545
Uniprot	Q6PIJ6
GeneID	81545;
Calculated MW	Refer to Figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200

## Background

This gene encodes a large protein that contains an F-box domain and may participate in protein ubiquitination. The encoded protein is a transcriptional co-activator of Krueppel-like factor 7 (Klf7). A heterozygous mutation in this gene was found in individuals with autosomal dominant distal hereditary motor neuropathy type IID. There is a pseudogene for this gene on chromosome 4. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only