

STRIP1 Polyclonal Antibody

Catalog No: #27382

Package Size: #27382-1 50ul #27382-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	STRIP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human STRIP1 (NP_149079.2).
Other Names	STRIP1;FAM40A;FAR11A
Accession No.	Uniprot:Q5VSL9GeneID:85369
Uniprot	Q5VSL9
GeneID	85369
Calculated MW	110kDa
SDS-PAGE MW	110kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:1000 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using STRIP1 antibody.

Background

This gene encodes a member of the striatin-interacting phosphatase and kinase complex, which is involved in localization of the Golgi body. The encoded protein participates in cytoskeletal organization. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only