

RBM4 Polyclonal Antibody

Catalog No: #27398

Package Size: #27398-1 50ul #27398-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RBM4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RBM4 (NP_002887.2).
Other Names	RBM4; LARK; RBM4A; ZCCHC21; ZCRB3A; RNA-binding protein 4
Accession No.	Swiss-Prot#:Q9BWF3NCBI Gene ID:5936
Uniprot	Q9BWF3
GeneID	5936;
Calculated MW	45kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RBM4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only