

RAB27B Polyclonal Antibody

Catalog No: #27414



Package Size: #27414-1 50ul #27414-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RAB27B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RAB27B (NP_004154.2).
Other Names	RAB27B; C25KG; ras-related protein Rab-27B
Accession No.	Swiss-Prot#:O00194NCBI Gene ID:5874
Uniprot	O00194
GeneID	5874;
Calculated MW	25kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RAB27B at 1:1000 dilution.

Background

Note: This product is for in vitro research use only