

TREML1 Polyclonal Antibody

Catalog No: #27425

Package Size: #27425-1 50ul #27425-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TREML1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TREML1 (NP_835468.1).
Other Names	TREML1; GLTL1825; PRO3438; TLT-1; TLT1; dJ238O23.3; trem-like transcript 1 protein
Accession No.	Swiss-Prot#:Q86YW5NCBI Gene ID:340205
Uniprot	Q86YW5
GeneID	340205;
Calculated MW	33kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TREML1 at 1:1000 dilution.

Background

This gene encodes a member of the triggering receptor expressed on myeloid cells-like (TREM) family. The encoded protein is a type 1 single Ig domain orphan receptor localized to the alpha-granule membranes of platelets. The encoded protein is involved in platelet aggregation, inflammation, and cellular activation and has been linked to Gray platelet syndrome. Alternative splicing results in multiple transcript variants

Note: This product is for in vitro research use only