

ATF7IP Polyclonal Antibody

Catalog No: #27442

Package Size: #27442-1 50ul #27442-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ATF7IP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ATF7IP (NP_060649.3).
Other Names	ATF7IP;AM;ATF-IP;MCAF;MCAF1;p621
Accession No.	Uniprot:Q6VMQ6GeneID:55729
Uniprot	Q6VMQ6
GeneID	55729
SDS-PAGE MW	250KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ATF7IP antibody.

Background

ATF7IP is a multifunctional nuclear protein that associates with heterochromatin. It can act as a transcriptional coactivator or corepressor depending upon its binding partners (summary by Liu et al., 2009 [PubMed 19106100]).[supplied by OMIM, Nov 2010]

Note: This product is for in vitro research use only