

CCDC131 Antibody

Catalog No: #45899

Package Size: #45899-1 50ul #45899-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CCDC131 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	CCDC131 Antibody detects endogenous levels of total CCDC131
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CCDC131
Target Name	CCDC131
Other Names	ZFC3H1, CCDC131, KIAA0546, PSRC2
Accession No.	Swiss-Prot#: O60293NCBI Gene ID:196441
Uniprot	O60293
GeneID	196441;
Calculated MW	226 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only