

RELL2 Polyclonal Antibody

Catalog No: #27487



Package Size: #27487-1 50ul #27487-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RELL2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human RELL2 (NP_776189.3).
Other Names	RELL2; C5orf16; RELT like 2
Accession No.	Swiss-Prot#:Q8NC24NCBI Gene ID:285613
Uniprot	Q8NC24
GeneID	285613;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IHC 1:50 - 1:200

Background

Note: This product is for in vitro research use only