

slc25a43 Polyclonal Antibody

Catalog No: #27494



Package Size: #27494-1 50ul #27494-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | slc25a43 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Human,Mouse,Rat |
| Immunogen Description | A synthetic peptide of human slc25a43 (NP_660348.2). |
| Other Names | SLC25A43; solute carrier family 25 member 43 |
| Accession No. | Swiss-Prot#:_NCBI Gene ID:194744 |
| GeneID | :194744 |
| Calculated MW | Refer to figures |
| Concentration | 1.9 mg/mL |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB □ 1:500 - 1:2000

Background

This gene encodes a member of the mitochondrial carrier family of proteins.

Note: This product is for in vitro research use only