

C1QL3 Antibody

Catalog No: #45903

Package Size: #45903-1 50ul #45903-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C1QL3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	C1QL3 Antibody detects endogenous levels of total C1QL3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human C1QL3
Target Name	C1QL3
Other Names	C1QL3, C1ql, CTRP13, Complement C1q-like protein 3, C1QTNF13, K100
Accession No.	Swiss-Prot#: Q5VWW1NCBI Gene ID:389941
Uniprot	Q5VWW1
GeneID	389941;
Calculated MW	27 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only