

## C1QL4 Antibody

Catalog No: #45904

Package Size: #45904-1 50ul #45904-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	C1QL4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	C1QL4 Antibody detects endogenous levels of total C1QL4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human C1QL4
Target Name	C1QL4
Other Names	C1QL4, CTRP11, C1QTNF11, Complement C1q-like protein 4
Accession No.	Swiss-Prot#: Q86Z23NCBI Gene ID:338761
Uniprot	Q86Z23
GeneID	338761;
Calculated MW	25 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only