

CNTNAP4 Antibody

Catalog No: #45908

Package Size: #45908-1 50ul #45908-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CNTNAP4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	CNTNAP4 Antibody detects endogenous levels of total CNTNAP4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CNTNAP4
Target Name	CNTNAP4
Other Names	CNTNAP4, CASPR4, KIAA1763
Accession No.	Swiss-Prot#: Q9C0A0NCBI Gene ID:85445
Uniprot	Q9C0A0
GeneID	85445;
Calculated MW	145 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only