

RFPL3 Polyclonal Antibody

Catalog No: #27678



Package Size: #27678-1 50ul #27678-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RFPL3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human RFPL3 (NP_006595.1).
Other Names	RFPL3; ret finger protein-like 3
Accession No.	Swiss-Prot#:O75679NCBI Gene ID:10738
Uniprot	O75679
GeneID	10738;
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RFPL3 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only