

USP25 Polyclonal Antibody

Catalog No: #27724



Package Size: #27724-1 50ul #27724-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	USP25 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human USP25 (NP_037528.3).
Other Names	USP25; USP21; ubiquitin specific peptidase 25
Accession No.	Swiss-Prot#:Q9UHP3NCBI Gene ID:29761
Uniprot	Q9UHP3
GeneID	29761;
Calculated MW	130kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using USP25 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only