

EVC Polyclonal Antibody

Catalog No: #27736



Package Size: #27736-1 50ul #27736-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EVC Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human EVC (NP_714928.1).
Other Names	EVC; DWF-1; EVC1; EVCL; EvC ciliary complex subunit 1
Accession No.	Swiss-Prot#:P57679NCBI Gene ID:2121
Uniprot	P57679
GeneID	2121;
Calculated MW	Refer to Figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:1000 - 1:2000

Background

This gene encodes a protein containing a leucine zipper and a transmembrane domain. This gene has been implicated in both Ellis-van Creveld syndrome (EvC) and Weyers acrodermal dysostosis.

Note: This product is for in vitro research use only