

## APOBEC3H Antibody

Catalog No: #45925

Package Size: #45925-1 50ul #45925-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	APOBEC3H Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	APOBEC3H Antibody detects endogenous levels of total APOBEC3H
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human APOBEC3H
Target Name	APOBEC3H
Other Names	APOBEC3H, APOBEC-related protein 10, ARP10, ARP-10
Accession No.	Swiss-Prot#: Q6NTF7NCBI Gene ID:164668
Uniprot	Q6NTF7
GeneID	164668;
Calculated MW	24 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only