

POLD4 Antibody

Catalog No: #45928

Package Size: #45928-1 50ul #45928-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	POLD4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	POLD4 Antibody detects endogenous levels of total POLD4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human POLD4
Target Name	POLD4
Other Names	POLDS
Accession No.	Swiss-Prot#: Q9HCU8NCBI Gene ID:57804
Uniprot	Q9HCU8
GeneID	57804;
Calculated MW	12 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only