

## DNA2L Antibody

Catalog No: #45931

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA2L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	DNA2L Antibody detects endogenous levels of total DNA2L
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human DNA2L
Target Name	DNA2L
Other Names	DNA2, DNA2-like helicase, DNA2L, HDNA2, KIAA0083
Accession No.	Swiss-Prot#: P51530NCBI Gene ID:1763
Uniprot	P51530
GeneID	1763;
Calculated MW	120 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only