

FAM62B Polyclonal Antibody

Catalog No: #27805



Package Size: #27805-1 50ul #27805-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM62B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human FAM62B (NP_065779.1).
Other Names	ESYT2; CHR2SYT; E-Syt2; FAM62B; extended synaptotagmin-2
Accession No.	Swiss-Prot#:A0FGR8NCBI Gene ID:57488
Uniprot	A0FGR8
GeneID	57488;
Calculated MW	102kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ESYT2 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only