

## PARVG Polyclonal Antibody

Catalog No: #27829

Package Size: #27829-1 50ul #27829-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PARVG Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human PARVG (NP_071424.1).
Other Names	PARVG; gamma-parvin
Accession No.	Swiss-Prot#:Q9HBI0NCBI Gene ID:64098
Uniprot	Q9HBI0
GeneID	64098;
Calculated MW	37kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of mouse testis, using PARVG at 1:3000 dilution.

## Background

---

Note: This product is for in vitro research use only