

RPAB4 Antibody

Catalog No: #45944

Package Size: #45944-1 50ul #45944-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RPAB4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	RPAB4 Antibody detects endogenous levels of total RPAB4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RPAB4
Target Name	RPAB4
Other Names	POLR2K, ABC10-alpha, HsRPB10a, HRPB7.0, RPB12, RPB7.0, RPABC4, RPB10alpha
Accession No.	Swiss-Prot#: P53803NCBI Gene ID:5440
Uniprot	P53803
GeneID	5440;
Calculated MW	7 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of POLR2K expression in DLD (A), SP2/O (B), PC12 (C) whole cell lysates.

Note: This product is for in vitro research use only