

IFIT3 Polyclonal Antibody

Catalog No: #27980



Package Size: #27980-1 50ul #27980-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	IFIT3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human IFIT3 (NP_001540.2).
Other Names	IFIT3; CIG-49; GARG-49; IFI60; IFIT4; IRG2; ISG60; P60; RIG-G; cig41; interferon-induced protein with tetratricopeptide repeats 3
Accession No.	Swiss-Prot#:O14879NCBI Gene ID:3437
Uniprot	O14879
GeneID	3437;
Calculated MW	60kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using IFIT3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only