

KIAA1456 Polyclonal Antibody

Catalog No: #28310



Package Size: #28310-1 50ul #28310-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	KIAA1456 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human KIAA1456 (NP_065895.2).
Other Names	KIAA1456; C8orf79; TRM9L
Accession No.	Swiss-Prot#:Q9P272NCBI Gene ID:57604
Uniprot	Q9P272
GeneID	57604;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IF□1:50 - 1:100

Background

Note: This product is for in vitro research use only