

## HECW1 Antibody

Catalog No: #45953



Package Size: #45953-1 50ul #45953-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HECW1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	HECW1 Antibody detects endogenous levels of total HECW1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HECW1
Target Name	HECW1
Other Names	KIAA0322, NEDL1
Accession No.	Swiss-Prot#: Q76N89NCBI Gene ID:23072
Uniprot	Q76N89
GeneID	23072;
Calculated MW	180kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only