

## MARCH1 Antibody

Catalog No: #45955

Package Size: #45955-1 50ul #45955-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description


Product Name	MARCH1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	MARCH1 Antibody detects endogenous levels of total 42430
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MARCH1
Target Name	42430
Other Names	MARCH1, MARCH-I, RING finger protein 171, RNF171
Accession No.	Swiss-Prot#: Q8TCQ1NCBI Gene ID:55016
Uniprot	Q8TCQ1
GeneID	55016;
Calculated MW	32 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images





Immunohistochemical analysis of paraffin-embedded human meningioma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

---

Note: This product is for in vitro research use only