

MARCH9 Antibody

Catalog No: #45959

Package Size: #45959-1 50ul #45959-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MARCH9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	MARCH9 Antibody detects endogenous levels of total 42438
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MARCH9
Target Name	42438
Other Names	MARCH9, RING finger protein 179, RNF179, MARCH-IX
Accession No.	Swiss-Prot#: Q86YJ5NCBI Gene ID:92979
Uniprot	Q86YJ5
GeneID	92979;
Calculated MW	38 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4o½ over night

Note: This product is for in vitro research use only