

RNF123 Antibody

Catalog No: #45962



Package Size: #45962-1 50ul #45962-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RNF123 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	RNF123 Antibody detects endogenous levels of total RNF123
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RNF123
Target Name	RNF123
Other Names	KPC1
Accession No.	Swiss-Prot#: Q5XPI4NCBI Gene ID:63891
Uniprot	Q5XPI4
GeneID	63891;
Calculated MW	149 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only