

RNF167 Antibody

Catalog No: #45963

Package Size: #45963-1 50ul #45963-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RNF167 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RNF167 Antibody detects endogenous levels of total RNF167
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RNF167
Target Name	RNF167
Other Names	RNF167, Ring finger protein 167, RING105, DKFZP566H073, LP2254
Accession No.	Swiss-Prot#: Q9H6Y7NCBI Gene ID:26001
Uniprot	Q9H6Y7
GeneID	26001;
Calculated MW	38 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4o½ over night

Note: This product is for in vitro research use only