

ZNF169 Polyclonal Antibody

Catalog No: #28460



Package Size: #28460-1 50ul #28460-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF169 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZNF169 (NP_919301.2).
Other Names	ZNF169; zinc finger protein 169
Accession No.	Swiss-Prot#:Q14929NCBI Gene ID:169841
Uniprot	Q14929
GeneID	169841;
Calculated MW	Refer to Figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF169 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only