

ENDOD1 Antibody

Catalog No: #45967

Package Size: #45967-1 50ul #45967-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ENDOD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ENDOD1 Antibody detects endogenous levels of total ENDOD1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ENDOD1
Target Name	ENDOD1
Other Names	ENDOD1, KIAA0830
Accession No.	Swiss-Prot#: O94919NCBI Gene ID:23052
Uniprot	O94919
GeneID	23052;
Calculated MW	55 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only