

NEIL2 Antibody

Catalog No: #45969

Package Size: #45969-1 50ul #45969-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NEIL2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	NEIL2 Antibody detects endogenous levels of total NEIL2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NEIL2
Target Name	NEIL2
Other Names	NEIL2, Nei homolog 2, NEH2, Nei like 2, Nei like 2 (E. coli), DNA glycosylase/AP lyase Neil2, Endonuclease 8-like 2, Endonuclease VIII-like 2, Nei-like protein 2, NEI2
Accession No.	Swiss-Prot#: Q969S2NCBI Gene ID:252969
Uniprot	Q969S2
GeneID	252969;
Calculated MW	37 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 40% over night

Note: This product is for in vitro research use only