

EDEM3 Antibody

Catalog No: #45972

Package Size: #45972-1 50ul #45972-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EDEM3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	EDEM3 Antibody detects endogenous levels of total EDEM3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human EDEM3
Target Name	EDEM3
Other Names	EDEM3, C1orf22
Accession No.	Swiss-Prot#: Q9BZQ6NCBI Gene ID:80267
Uniprot	Q9BZQ6
GeneID	80267;
Calculated MW	105kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of Hela lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only