

FANCM Antibody

Catalog No: #45977



Package Size: #45977-1 50ul #45977-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | FANCM Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Immunogen affinity purified |
| Applications | WB |
| Species Reactivity | Hu |
| Specificity | FANCM Antibody detects endogenous levels of total FANCM |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human FANCM |
| Target Name | FANCM |
| Other Names | FANCM, FAAP250, Fanconi anemia group M protein, KIAA1596, Protein FACM, Protein Hef ortholog |
| Accession No. | Swiss-Prot#: Q8IYD8NCBI Gene ID:57697 |
| Uniprot | Q8IYD8 |
| GeneID | 57697; |
| Calculated MW | 232 kDa |
| Concentration | 1.0mg/mL |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C |

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only