

FNBP4 Antibody

Catalog No: #45983

Package Size: #45983-1 50ul #45983-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FNBP4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	FNBP4 Antibody detects endogenous levels of total FNBP4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FNBP4
Target Name	FNBP4
Other Names	FBP30, KIAA1014
Accession No.	Swiss-Prot#: Q8N3X1NCBI Gene ID:23360
Uniprot	Q8N3X1
GeneID	23360;
Calculated MW	110 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only