

SLC19A1 Polyclonal Antibody

Catalog No: #28660



Package Size: #28660-1 50ul #28660-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC19A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SLC19A1 (NP_919231.1).
Other Names	SLC19A1; CHMD; FOLT; IFC1; REFC; RFC1; folate transporter 1
Accession No.	Swiss-Prot#:P41440NCBI Gene ID:6573
Uniprot	P41440
GeneID	6573;
Calculated MW	Predicted band size:65kDa
Concentration	0.2 mg/ml
Formulation	PBS with 0.01% thimerosal,50% glycerol,pH7.3.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Background

The membrane protein encoded by this gene is a transporter of folate and is involved in the regulation of intracellular concentrations of folate. Three transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only