

ALOX12 Polyclonal Antibody

Catalog No: #28678

Package Size: #28678-1 50ul #28678-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ALOX12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ALOX12 (NP_000688.2).
Other Names	ALOX12; 12-LOX; 12S-LOX; LOG12; arachidonate 12-lipoxygenase, 12S-type
Accession No.	Swiss-Prot#:P18054NCBI Gene ID:239
Uniprot	P18054
GeneID	239;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC□1:50 - 1:200IF□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded rat heart using ALOX12 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human colon carcinoma using ALOX12 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using ALOX12 at dilution of 1:100 (40x lens).



Immunofluorescence analysis of rat lung using ALOX12 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only