

GON4L Antibody

Catalog No: #45991

Package Size: #45991-1 50ul #45991-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GON4L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	GON4L Antibody detects endogenous levels of total GON4L
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GON4L
Target Name	GON4L
Other Names	GON4, KIAA1606
Accession No.	Swiss-Prot#: Q3T8J9NCBI Gene ID:54856
Uniprot	Q3T8J9
GeneID	54856;
Calculated MW	249 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of GON4L expression in HEK293T (A), K562 (B) whole cell lysates.

Note: This product is for in vitro research use only