

SRD5A1 Polyclonal Antibody

Catalog No: #28705

Package Size: #28705-1 50ul #28705-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SRD5A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Synthesized peptide derived from part region of human protein
Other Names	SRD5A1; S5AR 1; S5AR1; steroid 5 alpha-reductase 1
Accession No.	Swiss-Prot#:P18405NCBI Gene ID:6715
Uniprot	P18405
GeneID	6715;
Calculated MW	29kDa
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only