

## CAS9 Polyclonal Antibody

Catalog No: #28798

Package Size: #28798-1 50ul #28798-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CAS9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	WideRange
Immunogen Description	Recombinant protein of streptococcus pyogenes CAS9
Other Names	CAS9
Accession No.	Uniprot:/GeneID:901176
GeneID	901176
Calculated MW	150kDa
SDS-PAGE MW	170KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


## Application Details

WB □ 1:5000 - 1:10000 IF □ 1:50 - 1:200

## Images



Western blot analysis of extracts of 293T cells, using CAS9 antibody.



Immunofluorescence analysis of HeLa cells transfected with Cas9 (upper left) and untreated HeLa cells (upper right) using CAS9 Rabbit pAb (red, A14997) at dilution of 1:100. Blue: DAPI for nuclear staining.

---

Note: This product is for in vitro research use only