

## GNG3 Antibody

Catalog No: #46005



Package Size: #46005-1 50ul #46005-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GNG3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	GNG3 Antibody detects endogenous levels of total GNG3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GNG3
Target Name	GNG3
Other Names	GNG3, G gamma3, GNGT3, NBP gamma-3
Accession No.	Swiss-Prot#: P63215NCBI Gene ID:2785
Uniprot	P63215
GeneID	2785;
Calculated MW	8 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only