

HOXB3 Antibody

Catalog No: #46016

Package Size: #46016-1 50ul #46016-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HOXB3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	HOXB3 Antibody detects endogenous levels of total HOXB3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HOXB3
Target Name	HOXB3
Other Names	HOX2G
Accession No.	Swiss-Prot#: P14651NCBI Gene ID:3213
Uniprot	P14651
GeneID	3213;
Calculated MW	44 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only