

## HOXC13 Antibody

Catalog No: #46018

Package Size: #46018-1 50ul #46018-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HOXC13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	HOXC13 Antibody detects endogenous levels of total HOXC13
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HOXC13
Target Name	HOXC13
Other Names	HOXC13, Homeo box C13, Homeobox protein Hox-3G, Homeobox protein Hox-C13, Homeo box 3G, HOX3G, ECTD9, Homeobox C13, HOX3, NUP98/HOXC13
Accession No.	Swiss-Prot#: P31276NCBI Gene ID:3229
Uniprot	P31276
GeneID	3229;
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4o½ over night

---

Note: This product is for in vitro research use only