

MRGPRX4 antibody

Catalog No: #23165

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	MRGPRX4 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 258 and 322 of MRGPRX4
Target Name	MRGPRX4
Other Names	GPCR; MRGX4; SNSR6; MGC129753; MGC129754
Accession No.	Swiss-Prot:Q96LA9Gene ID:117196
Uniprot	Q96LA9
GeneID	117196;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 36kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: H1299
10% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Tw2.6 xenograft, using MRGPRX4 antibody at 1: 100 dilution.

Note: This product is for in vitro research use only